ANIOS SPS PREMIUM

ANIOS SPS Premium

Sommaire

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Regulation and general information

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*ANIOS SPS *Premium* is designed, produced and controlled by the Laboratoires ANIOS, certified by the AFAQ Organism under the number 1995/3723 in accordance with the ISO 9001 Quality System*

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*ANIOS SPS *Premium* is in compliance with the requirements of the European regulation related to biocides for disinfection of surfaces (Group 1 - PT 2)*

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*ANIOS SPS *Premium* formulation fulfils the requirements of Regulation (EC) No 648/2004 on detergents and its amendments*

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*ANIOS SPS *Premium* is labelled in accordance with the European regulation related to the classification and labelling of chemical products*

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*ANIOS SPS *Premium* meets our voluntary eco-design commitment included in our ANIOSAFE charter*

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Didecyldimethyl ammonium chloride

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Pure product ready to use :

- *Storage between +5°C and +35°C.*
- *Stability:*
 - *Not opened spray: 3 years from the manufacturing date indicated on the label.*
 - *Opened spray: preserving of the 3 years if the product is stored in a container closed by the sprayer.*

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Not concerned by biodegradability notion

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*Readily biodegradable substances and/or inherently biodegradable*³

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Substances which not fulfil the biodegradability criteria, or substances without any biodegradability data.

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ANIOS SPS PREMIUM contains at least 99.7 % of inorganic material and of biodegradable organic material.

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Environment

*This mixture does not present an environmental hazard. No known or foreseeable environmental damage under standard conditions of use.**

* Classification drawn up in accordance with European regulation related to the classification and labelling of chemical products.

PERSONAL PROTECTION UNDER THE CONDITIONS OF THE INTENDED USE

Refers exclusively to this product. If the product is handled and/or exposed simultaneously with other chemical agents, these must be taken into consideration when choosing personal safety equipment.

Respiratory protection

*In case of insufficient ventilation, with the risk of exceeding the TLVs/TWA**, wear suitable breathing apparatus.*

Particularly breathing apparatus, A2P2 type.

Hand protection

When handling this product, wear suitable gloves.

Nitrile, latex or vinyl gloves.

Gloves must be replaced immediately if they show signs of wear and tear.

Eye and face protection

Avoid contact with eyes.

Before handling, wear safety goggles with protective sides accordance with standard EN166.

Prescription glasses are not considered as protection.

Provide eye washes in workshops.

If necessary, the nearest water supply

** please refer to SDS §8 for TLVs (Threshold Limit Values)/TWA (Time-Weighted Average) of concerned substances present in the product

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The Safety Data Sheet (SDS) of this product is available on our Web site www.anios.com. In addition to visualization and/or impression, it is possible for you to download it.

The consultation of this SDS via internet enables you to receive later on by email any updated version of this document.

In the absence of an Internet access, the SDS remains available on request near our laboratories (tel.: +33 3 20 67 67 67 - fax: +33 3 20 67 67 68)

* an access is also possible at the site <http://www.quickfds.com>