ANIOSYME X3

ANIOSYME X3

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Regulation and general information

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ANIOSYME X3 is designed, produced and controlled by the Laboratoires ANIOS, certified by the AFAQ Organism under the number 1995/3723 in accordance with the ISO 9001 Quality System

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ANIOSYME X3 has the CE mark for medical devices in the category class IIb, in accordance with the Directive 93/42/CE

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ANIOSYME X3 formulation fulfils the requirements of Regulation (EC) No 648/2004 on detergents and its amendments

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ANIOSYME X3 is labelled in accordance with the European regulation related to the classification and labelling of chemical products

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ANIOSYME X3 meets our voluntary eco-design commitment included in our ANIOSAFE charter

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Digluconate de chlorhexidine 0.30%
Chlorhexidine digluconate

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Other ingredients

Complexe enzymatique (protéase, lipase et amylase)
Enzymatic complex (protease, lipase and amylase)

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Pure product ready to use :

- Storage between +5°C and +35°C
- Stability :
 - Not opened bottle : 24 months from the manufacturing date indicated on the label.
 - Opened bottle : maintenance of the 24 months in condition if closing the bottle after use.

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Diluted product

Dilution: 0.5% - 5mL per liter of water

Conservation of the diluted solution in a closed container (cap or spray gun) : 7 days.

Conservation of the diluted solution in an opened container : 24 hours (without use).

Any diluted preparation will be eliminated after use.

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ANIOSYME X3 antimicrobial efficacy is tested against monobacterial biofilm with Pseudomonas aeruginosa, artificially developed inside a Tygon® tube, according to the ISO/TC 15883 – 5 : 2006 standard.

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After the passing of the 0.5% ANIOSYME X3 diluted solution during 5 minutes, followed by rinsing with distilled water during 2 minutes, number of viable bacteria fixed by cm² of carrier and concentration of proteins are determined inside the Tygon® tube.

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Concentration in proteins determined in biofilm

	Concentration en protéines $\mu\text{g}/\text{cm}^2$ <i>Concentration in proteins $\mu\text{g}/\text{cm}^2$</i>	% de réduction <i>Reduction %</i>
Tube témoin <i>Control tube</i>	17,4	
ANIOSYME X3 (0.5%)	2,8	<u>84,2%</u>

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
Present surfactants distinguish themselves by hydrophilic and hydrophobic properties, which ensures performances against proteic and glucidic soilings on one hand, and greasy soiling (lipidic) on another hand, by solubilisation and/or emulsification.

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The tri-enzyme pre-disinfectant detergent ANIOSYME X3, diluted at 0.5 % in tap water, presents a High Degreasing Cleaning Power against a blood based soiling fixed on carriers made of brass (hydrophilic character).

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Proteases, lipases and amylases act in the same way as surfactants, degrading, by hydrolysis, organic compounds not much or not soluble in soluble components with lowest molecular weight.

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The enzymatic activity of ANIOSYME X3 solution is demonstrated after accelerated ageing at +30°C during 56 days, then diluted at 0.5% and putting into contact with the substrate during 15 minutes at room temperature.

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
The realization of an electrophoresis on polyacrylamide gel makes it possible to highlight the proteasic activity.

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Lipasic and amylasic activities of ANIOSYME X3 are demonstrated identically by thin layer chromatography method.

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
Enzymatic stability of diluted solution (1)

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Enzymatic stability of diluted solution (2)

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The proteasic activity of a 0.5% ANIOSYME X3 solution is demonstrated after storage up to 72 hours at room temperature, by putting it into contact during 15 minutes with the human albumin substrate.

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
It is demonstrated, by electrochemical study of pitting corrosion according to the NF S 94-402-1 standard (May 2004), that ANIOSYME X3 does not present pitting corrosive character with stainless steel Z30 Cr13 in the conditions of use recommended by Laboratoires ANIOS.

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Not concerned by biodegradability notion

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*Readily biodegradable substances and/or inherently biodegradable*³

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Substances which not fulfil the biodegradability criteria, or substances without any biodegradability data.

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ANIOSYME X3 contains at least 99.8 % of inorganic material and of biodegradable organic material.

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HAZARDS IDENTIFICATION AT DILUTION OF USE*

Health

Not classified at dilution of use in terms of health hazards.

Environment

Not classified at dilution of use in terms of environment hazards.

* Classification drawn up in accordance with European regulation related to the classification and labelling of chemical products.

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PERSONAL PROTECTION AT THE APPLICATION RATE UNDER THE CONDITIONS OF THE INTENDED USE

Refers exclusively to this product. If the product is handled and/or exposed simultaneously with other chemical agents, these must be taken into consideration when choosing personal safety equipment.

Respiratory protection

In case of insufficient ventilation, with the risk of exceeding the TLVs (Threshold Limit values)/TWA, wear suitable breathing apparatus (mask which filters organic vapours - type A protection).

Hand protection

When handling, wear suitable gloves (nitrile, latex or vinyle).

Eye and face protection

*Avoid all contact with eyes.
Water source nearby.*

** please refer to SDS §8 for TWA (Time-Weighted Average)/STEL (Short-term exposure limit) of concerned substances present in the product

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The Safety Data Sheet (SDS) of this product is available on our Web site www.anios.com. In addition to visualization and/or impression, it is possible for you to download it.

The consultation of this SDS via internet enables you to receive later on by email any updated version of this document.

In the absence of an Internet access, the SDS remains available on request near our laboratories (tel.: +33 3 20 67 67 67 - fax: +33 3 20 67 67 68)

* an access is also possible at the site <http://www.quickfds.com>