Quickie Life

Directions for use
Notice d'utilisation
Manuale d'uso
Gebruikershandleiding
Manual de usuario
Instruções para Utilização
Brukerveiledning
Bruksanvisning
Brugervejledning

QUICKIE®

**INNOVATE.™
BUILD FOR THE FUTURE**

For visually impaired people, you will find information on our website
www.SunriseMedical.com.

Foreword

Dear Customer,

We are very pleased that you have decided in favour of a high-quality SUNRISE MEDICAL product.

This Owner's manual will provide numerous tips and ideas so that your new wheelchair can become a trustworthy and reliable partner in your life.

Maintaining close ties to our customers means a lot to us at Sunrise Medical. We like to keep you up-to-date about new and current developments at our company. Keeping close to our customers means: fast service, as little red tape as possible, working closely with customers. When you need replacement parts or accessories, or if you just have a question about your wheelchair – we are there for you.

We want you to be satisfied with our products and service. At Sunrise Medical we are constantly working to develop our products further. For this reason, changes can occur in our palette of products with regard to form, technology, and equipment. Consequently, no claims can be construed from the data or pictures contained in this owner's manual.

The management system of SUNRISE MEDICAL is certified to DIN EN ISO 9001, ISO 13485 and ISO 14001.

As the manufacturer, SUNRISE MEDICAL, declares that the lightweight wheelchairs conform to the 93/42/EWG 47/70/2010EWG guideline.

In addition, they fulfil the performance requirements for the "Crash Test" to ISO 7176-19

Please contact your local, authorised SUNRISE MEDICAL dealer if you have any questions regarding the use, maintenance or safety of your wheelchair.

In case there is no authorised dealer in your area, or you have any questions, contact Sunrise Medical either in writing or by telephone.

Sunrise Medical Limited
High Street, Wollaston
West Midlands DY8 4PS
Tel.: +44 (0) 1384446666
Fax: +44 (0) 1384446644
E-mail: enquiries@sunrisemedical.co.uk
Internet: www.sunrisemedical.co.uk

IMPORTANT:
DO NOT USE YOUR WHEELCHAIR UNTIL THIS MANUAL HAS BEEN READ AND UNDERSTOOD.

Table of contents

Foreword for Wheelchairs

Safety Tips for Wheelchairs

Overview	5
Guarantee	5
Handling for transportation	6-7

Options	
Handling	9
Step Tubes	9
Brakes	9-10
Footboard	10
Footrest	10
Castor	11
Perforated plate	11
Camber	11
Backrest	11
Sideguard	12
Push handle	12
Anti-tip tubes	12
Crutch holder	12
Safety belt	13
Stabilising bar	14
Transit wheels	14

Tyres and Mounting 14

Trouble Shooting 14

Disposal / Recycling of materials 14

Maintenance and care 15

Nameplate 15

Technical Specifications 15-16

Torque 17

Use

Wheelchairs are exclusively for a user who is unable to walk or has limited mobility, for their own personal use, when self-propelling and being moved by a third party (pushed by attendant) within the home and outdoors.

The maximum weight limit (includes both the user and any weight of accessories fitted to the wheelchair) is marked on the serial number label, which is affixed to the cross-brace or the stabilising bar under the seat.

Warranty can only be issued if the product is used under the specified conditions and for the intended purpose.

The expected life of the wheelchair is 5 years.

NO NON-STANDARD PARTS are to be fitted, other than those officially approved by Sunrise Medical.

Use

The variety of fitting variants as well as the modular design means that it can be used by those who cannot walk or have limited mobility because of:

- Paralysis
- Loss of extremity (leg amputation)
- Extremity defect/deformity
- Joint contractures/joint injuries
- Illnesses such as heart and circulation deficiencies,
- disturbance of equilibrium or cachexia as well as for elderly people who still have strength in the upper body.

When considering provision, please also note the body size, weight, physical and psychological constitution, the age of the person, living conditions and environment.

General safety notes and driving limits

The construction of this wheelchair, as well as the technology used, are designed to give the maximum level of safety. International safety standards currently in force have either been fulfilled or exceeded in parts. Nevertheless, users may put themselves at risk by improperly using their wheelchairs. For your own safety, the following rules must absolutely be observed.

Unprofessional or erroneous changes or adjustments increase the risk of accident. As a wheelchair user, you are also part of the daily traffic on streets and pavements, just like anyone else. We would like to remind you that you are, therefore, subject to any and all traffic laws.

Be careful during your first ride in this wheelchair. Get to know your wheelchair.

Before each use, the following should be checked:

- Quick-release axles on the rear wheels
- Velcro on seats and backrests
- Tyres, tyre pressure, and wheel locks.

Before changing any of the adjustments of this wheelchair, it is important to read the corresponding section of the Owner's manual.

Maximum user weight: 125 kg. Please note the weight information for lighter weight options, which are quoted separately.

This wheelchair should be used exclusively to convey one person at a time. Any other use does not comply with the intended purpose.

Only permanently installed seats and safety belts will offer sufficient protection in hazardous situations.

When vehicles are in motion, unoccupied wheelchairs should be secured using appropriate means.

Before changing any of the adjustments of this wheelchair, it is important to read the corresponding section of the Owner's manual.

It is possible that potholes or uneven ground could cause this wheelchair to tip over, especially when riding uphill or downhill.

When riding over a step or up an incline frontally, the body should be leaning forward.

Danger!

Secure your wheelchair on uneven ground or when transferring (e. g. into a car) by using the brakes.

Danger!

For thigh amputees you must position the rear wheels towards the back and/or use anti-tip tubes.

Danger!

With active centre of gravity adjustment, the wheelchair may tip over, so anti-tip tubes should be used.

Danger!

Before setting off, check that your tyre pressure is correct. The correct tyre pressure is printed on the edge of the wheel, for rear wheels it should be at least 3.5 bar (350 kPa).

All brakes will only work if there is sufficient tyre pressure and if the correct setting has been made (see the Chapter on "Brakes").

Danger!

The wheel locks are not intended to brake your wheelchair whilst it is moving. They are only there to ensure that your wheelchair does not begin rolling away. When you stop on uneven ground/floors, you should always use the wheel locks to prevent such rolling. Always apply both wheel locks; otherwise, your wheelchair could tip over. Always make use of elevators and ramps. If these are not available, you should request the help of an attendant. They should only grip the wheelchair at securely mounted parts. If anti-tip tubes are fitted these must be folded away. A wheelchair should never be lifted with an occupant; it should only be pushed.

Danger!

Do not lift or carry the wheelchair by the backtubes or the pushhandles.

Danger!

If the seat and back sling are damaged, you must replace them immediately.

Danger!

Be careful with fire, in particular with lighted cigarettes, seat and back slings can catch fire. seat and back slings can catch fire.

Danger!

To avoid injuries to the hand when operating the wheelchair, grip between the rear wheel and the brake.

Danger!


If a journey in a specially fitted vehicle for disabled people is necessary, we recommend that the user should use the seats in the vehicle and the appropriate restraint systems.

Danger!

Always make sure that the quick-release axles on the rear wheels are set properly. If the button on the quick-release axle is not pressed in, the rear wheel cannot be removed.

Danger!

In particular when using lightweight metal handdrims, fingers will

easily become hot when braking from a high speed or on long inclines.

Danger!

If the wheelchair is subject to direct sunlight for a long period of time, then parts of the wheelchair (e.g. frame, legrests, brakes and sideguard) may become hot (>41°C).

Note!

When using the chair outside, always wear leather gloves which will increase your grip and protect the fingers from dirt and overheating.

Note!

The effect of the brake as well as the general driving characteristics are dependent on tyre pressure. The wheelchair is significantly lighter and easier to manoeuvre when the rear wheels are pumped up correctly and both wheels have the same pressure.

Note!

Make sure that your tyres have sufficient tread!

Please note that you are subject to any and all traffic laws when driving in public traffic.

Note!

Before using the wheelchair for the first time, the user must receive the necessary instruction.

Note!

When it is dark, please wear light clothing or clothing with reflectors, so that you can be seen more easily. Make sure that the reflectors on the side and back of the wheelchair are clearly visible. We would also recommend that you fit an active light.

Note!

Do not use the wheelchair on slopes > 10°.

Do not use the wheelchair on muddy or icy ground.

Only use the wheelchair in areas which are permitted for pedestrians.

Note!

When working with or making adjustments to your wheelchair, be careful not to get your fingers caught!

Please keep a note of your local service agent's address and telephone number in the space provided.

In the event of a breakdown, contact them and try to give all relevant details so they can help you quickly.

The wheelchairs shown and described in this user guide may not correspond in every detail exactly to your own model. However, all instructions are completely relevant, regardless of possible detail differences.

The manufacturer reserves the right to alter without notice any weights, measurements or other technical data shown in this manual. All figures, measurements and capacities shown in this manual are approximate, and do not constitute specifications.

Operating conditions:

The wheelchair can be used indoors as well as outdoors on dry, relatively smooth surfaces. It should not be used in heavy rain, on snow-covered or slippery surfaces.

Guarantee

Guarantee

THIS DOES NOT AFFECT YOUR LEGAL RIGHTS IN ANY WAY.

Guarantee conditions

1) Repair or replacement is carried out by the authorised Sunrise Medical dealer.

2) To fulfil the guarantee conditions, should servicing need to be carried out on your wheelchair under this agreement, contact the designated Sunrise Medical customer service agent immediately, with precise details on the type of difficulty. Should you be using the wheelchair outside the area covered by the designated Sunrise Medical customer service agent, the work will be carried out under "guarantee conditions" by another agency as designated by the manufacturer.

3) Should a part or parts of the wheelchair require repair or replacement within 24 months (5 years for frame and cross-brace) after transfer of ownership to the original purchaser, and provided that this person is still the owner of the wheelchair, as a result of a specific manufacturing or material defect, the part or parts will be repaired or replaced free of charge if the wheelchair is returned to the authorised customer service agent.

Note: This guarantee cannot be transferred.

4) The guarantee also covers all repaired or replaced parts for the remaining period of the guarantee for the wheelchair.

5) For spare parts which are fitted after the start of the original guarantee, we give a further 24-month guarantee.

6) Consumable parts are normally excluded from the guarantee, except in the case that premature wear of the part is the direct result of a manufacturing fault. These parts include, amongst others, upholstery, tyres, inner tubes and similar parts.

7) The guarantee conditions above cover all product parts for models which were purchased at full sales price.

8) Normally we do not accept responsibility if a repair or replacement of the wheelchair is required for one of the following reasons:

a) The product or part has not been maintained or serviced in accordance with the manufacturer's recommendations as shown in the User Instructions and/or the Service Instructions. Accessories have been used which are not specified as original accessories.

b) The wheelchair or a part of the wheelchair was damaged through neglect, accident or improper use.

c) Alterations to the wheelchair or parts, which are not in accordance with the manufacturer's specifications or the carrying out of repairs before informing the customer service agent.

Transport

Transportation of your wheelchair within a vehicle

A wheelchair secured in a vehicle will not provide the equivalent of safety and security of a vehicle seating system. It is always recommended that the wheelchair user is transferred to a seat once in the vehicle. However, it is recognised, that it is not always possible for the wheelchair user to be transferred. In such cases, for example, when the wheelchair user must be transported while sitting in the wheelchair, then the following procedure should be followed:

1. Make sure that the vehicle is correctly equipped for the transportation of wheelchair and user, and that it is possible to get the wheelchair in and out of the vehicle. The vehicle should have the appropriate floor strength to take the combined weight of the wheelchair user, wheelchair and accessories.
2. There should be sufficient space available around the wheelchair, so that the restraint belt and fasteners are accessible and can be tightened or released again.
3. The occupied wheelchair must be located in a forward facing position and be secured by the wheelchair tie down system, and the wheelchair user must be secured by the belt system (restraint system according to WTORS in accordance with ISO 10542 or SAE J2249) and in accordance with the manufacturer's instructions (WTORS).
4. The wheelchair's use in other positions within a vehicle has not been tested e.g. transportation in a side facing position must not be carried out under any circumstances (Fig. a).

There is a risk of serious injury or death if this advice is ignored!



5. The wheelchair should be secured by a tie down restraint system according to ISO 10542 or SAE J2249 with non-adjustable front

straps and adjustable rear straps using Karabiner clips / S hooks and buckle fittings. This tie down restraint system consists of four individual belts, which are fixed to the four corners of the wheelchair.

6. The tie down restraint system should be fixed to the main frame of the wheelchair (see the diagram on the following page), under no circumstances should it be fixed to wheelchair components or accessory parts, i.e. not to the spokes of the wheel, the brakes or the legrests.
7. The tie down restraints should be attached as close as possible at an angle of 45 degrees and tightened securely in accordance with the manufacturer's instructions.

8. No changes should be made to the fixing points on the wheelchair or the structural parts, the frame or other wheelchair components without prior agreement from the manufacturer. If this advice is not followed, then a mobility aid from Sunrise Medical may not be transported in a vehicle.

9. To secure the wheelchair occupant and to reduce the risk of possible head and chest injury from impact with vehicle components and potential consequential injuries for the wheelchair user and other vehicle occupants, pelvic and upper torso belts must be used. (Fig. b) The upper torso belt must be fixed to the 'B' pillar of the vehicle – if this advice is not closely followed it could lead to an increased risk of abdominal injuries to the wheelchair user.

10. A head restraint suitable for transportation (see label on headrest) must be available and fitted accordingly during transportation.

11. Belts to support posture (pelvic restraint belts etc.) may only be used to secure the wheelchair user in a moving vehicle, if they meet the requirements of ISO 7176-19:2001 or SAE J2249 (see label).

12. The safety of the wheelchair user during transportation depends on the care and attention of the person who fits the restraint belts. The person who undertakes this type of task should have been trained in the use of such restraint systems.

13. Where possible, all additional accessories should be removed and stored safely, i.e. not on the wheelchair.

This includes for example:

- crutches, loose seat cushions, tray tables.
14. Swingaway or height-adjustable legrests should not be used in the raised position, if the wheelchair and wheelchair user is being transported and the wheelchair is clamped by the tie down restraint system and the wheelchair user is secured by seat belts.


15. Angle-adjustable backrests should be moved to the upright position.

16. The manual brakes must always be firmly applied.

17. Safety belts must be fixed to the "B" pillar of the vehicle and must not be held away from the body by wheelchair components or parts such as armrests or wheels.

Safety belts for wheelchair user - Instructions


1. The pelvic restraint belt must be fitted low across the lap, so that the angle of the pelvic belt is within the preferred zone of 30 - 75 degrees to the horizontal.


A steeper (greater) angle within the preferred zone is desirable i.e. closer to, but never exceeding 75°. (Fig. c)

2. The upper torso belt must be routed over the shoulder and diagonally across the upper torso - see Fig. d and e. Restraint belts must be adjusted as tightly as possible, whilst still being comfortable. The belts themselves must not be twisted when in use.

The upper torso belt must be routed over the shoulder and diagonally across the upper torso - see Fig. d and e.

3. The securement points on the wheelchair are the front inner side of the frame directly above the castors and on the rear of the frame. The belts are routed round the side frame at the intersection between the horizontal and vertical frame supports. (See Fig. 24-28)


4. The restraint belt symbol (Fig. f) on the wheelchair frame shows the position of the wheelchair restraint belt.
Once the front belts have been fitted, the belts are tightened to secure the wheelchair.

Transportability – Positioning of the wheelchair tie down belts on the mobility product

Wheelchair secured with front and rear wheelchair tie down belts (Fig. g).

The positioning of the restraint belts is shown in detail below.
Position of the front wheelchair tie down and the tie down label for the fixed frame and folding wheelchairs (Fig. h and i).


Position of the rear wheelchair tie down and tie down label on the wheelchair (Fig. j).

Wheelchair components

Wheelchairs:

1. Push handles
2. Back sling
3. Sideguard
4. Seat sling
5. Footrest
6. Castors
7. Footboard
8. Fork
9. Quick-release axle
10. Wheel locks
11. Handrim
12. Rear wheel

Handling

Folding up

First remove the seat cushion from the wheelchair and flip-up the footboard (continuous) or footboards.

Take hold of the sling or the seat tubes in the middle, from below and pull upwards. The chair will fold up.

To fold your wheelchair so that it is as small as possible, e.g., to put into a car, you can remove the footrests (depending on the model). To do this, push the release lever from the outside and swing the footrest to the side, then pull the footrest out of the frame tube.

Fig. 1

Fig. 2

Unfolding

Before unfolding, undo the folding bracket

Press down on the seat tube (see figure). The wheelchair will unfold.

Then press the seat tube into the support. This process is easier if you tip the chair a little, so that there is no load on one of the rear wheels and then the chair will unfold more easily.

Caution!

Risk of trapping fingers.

Quick-release axles on rear wheel

The rear wheels are equipped with quick-release axles. The wheels can, therefore, be installed or removed without using tools. To remove a wheel, simply depress the quick-release button on the axle (1) and pull it out.

Fig. 3

CAUTION!

Hold the quick-release button on the axle depressed when inserting the axle into the frame to mount the rear wheels. Release the button to lock the wheel in place. The quick-release button should snap back to its original position.

Options - Step Tubes

Step Tubes

Step tubes are used by attendants to push a wheelchair over an obstacle. Simply step on the tube to push a wheelchair, for example, over a kerb or step.

Fig. 4

Options - Brakes

Wheel locks


Your wheelchair is equipped with two wheel locks. They are applied directly against the tyres. To engage, press both wheel lock levers forward against the stops (Fig. 5).

To release the wheel lock, pull the levers back to their original positions.

Braking power will decrease with:

- Worn tyre tread
- Tyre pressure that is too low
- Wet tyres
- Improperly adjusted wheel locks

Fig. 5

The wheel locks have not been designed to be used as brakes for a moving wheelchair. The wheel locks should therefore never be used to brake a moving wheelchair. Always use the handrims for braking. Make sure that the interval between the tyres and wheel locks complies with the given specifications. To re-adjust, loosen screw (1) and set the appropriate interval. Tighten screw (see the page on torque).

CAUTION!

After each adjustment of the rear wheels, check the interval to the wheel locks and re-adjust if necessary.

Options - Brakes

Brake lever extension

The longer lever helps to minimise the effort needed to set the wheel locks. The brake lever extension is screwed to the brakes. By raising this, it can be flipped forward.

CAUTION!

Incorrect mounting of the wheel lock will result in a higher effort to operate. This might cause the wheel lock extension lever to break.

Leaning onto the wheel lock extension lever while transferring will cause the lever to break. There is a risk it could break. Moving parts may be made dirty by spray water.

Options - Footboard

Footrest - Plastic

Footplates:

The footplates can be flipped up to make it easier to transfer to/from the chair.

Lower leg length:

For users with longer legs, the leg support can be lowered to the appropriate level by using the screws. Loosen the screws, push the tubes to the desired position and then re-tighten the screws (see section on torque).

There must always be a minimum gap of 2.5 cm between the footplates and the ground. (Fig. 11)

Angle-adjustable footboard (7.22)

It can be adjusted to alter the angle to the ground. Release screw, pull it inwards, set the desired angle and then push it on. After making the adjustment, re-tighten the screw. Make sure that after any adjustment work, all screws are tightened correctly (see the page on torque).

Ensure that you maintain the minimum gap to the floor (2.5 cm) (Fig. 12).

Footrests and latch

The legrests can be swung inwards under the seat sling or outwards.

When fitting the legrests, the footrests must be pointing inwards or outwards. Then rotate the legrest inwards, until it locks into place. To remove, pull the lever (1), swing the footplate inwards or outwards and then lift the legrest out. Make sure that the legrest is locked into place correctly (Fig. 13).

CAUTION:

The legrests should not be used to lift or to carry the wheelchair.

CAUTION:

When getting on or off the wheelchair, do not use the footboards. These should be flipped up beforehand and be swung to the outside as far as possible.

Flip-up footrest

To elevate:

Take the load off the footrest and release the hand wheel. Pull the footrest upwards and set the desired height. Secure the desired height by tightening the hand wheel.

To lower:

Take the load off the footrest by lifting the lower leg and release the interlock with the hand wheel (Fig. 14).

In this way you can swing the footrest downwards. As soon as you have reached the desired position, secure the interlock with the hand wheel. You can set the lower leg length by releasing the screws, this is infinitely adjustable. (Fig. 15).

CAUTION:

Keep hands clear of the adjustment mechanism between the frame and the movable parts of the footrest while elevating or lowering the footrest.

CAUTION:

The footrests should not be used to lift or to carry the wheelchair.

CAUTION: Risk of trapping fingers!

When moving the footrest up or down, do not put your fingers in the adjusting mechanism between the moving parts of the footrest.

Options - Footboard

Footrest - Aluminium

There are several footboards available.

These can be flipped up to make it easier to transfer to/from the chair.

Lower leg length

By undoing the clamp screw (1), you can adjust the lower leg length. Loosen the clamp screw, adjust the tube with the footboard to the desired position.

Tighten clamp screws again (see the page on torque). Ensure that the plastic push joint is located under the clamp screw in the correct position.


Ensure that the plastic push joint is located under the clamp screw in the correct position.

Angle-adjustable footboard The footboard can be flipped up to make it easier to transfer to/from the chair.

It can be adjusted to alter the angle to the ground.

Tighten screw (1) on the outside. By loosening the clips (2) the footboard can be moved to three positions both forwards and backwards. By adjusting the set screw (3) the footboard can be changed to a horizontal position. To do this, the footboard must be flipped up. Make sure that after any adjustment work, all screws are tightened correctly (see the page on torque). Ensure that you maintain the minimum gap to the floor (2.5 cm).

Options - Castor


Castor, Castor adapter, Castor fork

From time to time the wheelchair may veer slightly to the right or left, or the castors may flutter. This may be caused by the following:

- The forward or reverse wheel motion has not been set properly.
- The castor angle has not been adjusted properly.
- Castor and/or rear wheel air pressure is incorrect; wheels do not turn smoothly.

The optimum adjustment of the castors is required so that the wheelchair runs in a straight line. The castors should always be adjusted by an authorised dealer.

The castor plates must be re-adjusted, and the wheel locks must be checked any time the rear wheel position has been altered. Check the setting of the castor adapter (1) by placing a set-square next to it and looking at it from the front. Make sure that the adapter is at right angles to the floor. By undoing the screws (2) and moving the adapter, you can then re-adjust the wheels to their optimum position, by adjusting the corresponding teeth segments. Make sure that after any adjustment work, all screws are tightened correctly (see the page on torque).

Options - Camber


Camber

The camber of the chair can be adjusted from 1° to 4°.

Doing this gives increased stability against the wheelchair falling to the side, and will also improve the wheelchair's manoeuvrability.


Loosen the screws on the perforated plate (1). Then adjust the adaptor (2) to the desired angle setting. Then tighten the screws again (see the page on torque).

Setting the wheel camber increases the overall width of the wheelchair (each degree adds on approx. 1 cm).

CAUTION!

If necessary, the brakes and castors must be re-adjusted.


Options - Backrest

Angle-adjustable fixed back


The angle-adjustable backrest can be set to 7 positions (9° forwards, 0°, 3°, 6° and 9° backwards). By loosening and removing the screw (1) the angle can be adjusted in stages of 3°.

For your own safety, adjustment work should always be carried out by your dealer.

Adjustable back sling

The adjustable back sling can be adjusted for tension by using several straps. The back sling upholstery can be accessed from the inside via an opening and can be padded to suit individual tastes.


Options - Perforated plate

Perforated plate

The position of the centre of gravity in relation to the rear wheel is an important factor for easy and comfortable driving. Several positions are possible. If you move the perforated plate to the rear (1) in the frame holes (2), the chair is in a safe position, but at the same time is not so manoeuvrable, as it is when the perforated plate is positioned further forward.

By moving the axle adapter (3) in the perforated plate, the seat height of the chair can be adjusted. Make sure that after any adjustment work, all screws are tightened correctly (see the page on torque).

CAUTION!

If necessary, the brakes and castors must be re-adjusted.

Options - Sideguard

Standard sideguard, flip-up, removable with short or long armrests


Sideguard, flip-up, removable with short or long armrests, height-adjustable The sideguard with rounded off front edges will allow you to get up close to a table.

The armrest can be height-adjusted as follows.

Push the lever upwards and move the armrest (2) to the desired height. Release the lever and press the armrest (2) downwards until you hear it click into place. Always check that the sideguards are correctly engaged. To flip the armrest up, pull the lever (1) so that the sideguard is released.

CAUTION!

Neither the sideguards nor the armrests are to be used for lifting or carrying the wheelchair.

Sideguard with clothes protector

The clothes protector prevents clothes getting dirty from spray water. You can set the position in relation to the rear wheel by moving the sideguard.

To do this, remove the screws (1 and 2). After setting to the desired position, re-tighten the screws (see the page on torque).

CAUTION!

Neither the sideguards nor the armrests are to be used for lifting or carrying the wheelchair

Options - Push handle

Height-adjustable push handles


These handles are secured with pins to prevent them from sliding out unintentionally. Opening the quick-release lever makes it possible to adjust the height of the push handles to meet your individual needs. As you move the lever, you will hear a locking mechanism; you may now easily position the push handle as desired. The nut on the tension lever determines how tightly the push handles are clamped into place. If the nut is loose after adjusting the tension lever, the push handle will also be too loose. Turn the push handle from side to side before use to make sure that it is clamped securely enough into place. After adjusting the handle height, always clamp the tension lever securely into place. If the lever is not secure, injuries could result when ascending stairs.

Options - Push handle

Fold-down push handles


If the push handles are not in use, they can be folded down depressing the button. When they are needed again, simply flip them back up until they click into place.

Options - Anti-tip tubes

Anti-tip tubes

The anti-tip tubes (1) provide additional safety for inexperienced users when they are still learning how to operate their wheelchairs. Anti-tip tubes (1) prevent a wheelchair from tipping over backwards. The anti-tip tube (1) can be swung downwards by pushing it or rotated forwards. There should be a gap of 3 to 5 cm to the floor. To go up or down a step (e.g. kerb) the anti-tip tubes must be swung forwards, to prevent them from touching the ground.

Options - Crutch holder

Crutch holder

This device permits crutches to be transported directly on the wheelchair. It has a Velcro loop to fasten crutches or other aids.

CAUTION!

Never try to use or even remove the crutches or other aids while moving.

Options – Pelvic restraint belt

Before using the wheelchair make sure that the pelvic restraint belt is fitted.

Fig. 28

The pelvic restraint belt must be checked on a daily basis to ensure it is free from any obstruction or adverse wear.

Always make sure that the pelvic restraint belt is correctly secured and adjusted prior to use. If the strap is too loose, this could cause the user of the wheelchair to slip down and risk suffocation or serious injury.

Fig. 29

The pelvic restraint belt is fitted to the wheelchair as shown in the illustrations. The seat belt comprises two halves. They are fitted using the existing seat stay retaining bolt fitted through the eyelet on the belt. The belt is routed under the rear of the side panel. (Fig. 28)

Adjust the belt position so buckles are in the centre of the seat. (Fig. 29)


Adjust the pelvic restraint belt to suit the user's needs as follows:

To reduce the belt length	To increase the belt length
	
Feed free belt back through male buckle and slide adjusters. Ensure belt is not looped at male buckle.	Feed free belt through slide adjusters and male buckle to provide more belt length. Fig. 30 

When fastened check space between the pelvic restraint belt and user. When correctly adjusted it should not be possible to insert more than the flat of the hand between the pelvic restraint belt and the user. (Fig. 30)

The pelvic restraint belt should be fixed so that the belt sits at an angle of 45 degrees across the user's pelvis. The user should be upright and be as far back as possible in the seat when correctly adjusted. The pelvic restraint belt should not allow the user to slip down in the seat. (Fig. 31)

Fig. 31

To fasten buckle:

Firmly push male buckle into female buckle.

To release belt:

Press exposed sides of male buckle together at the same time and push to the middle, while at the same time pull them gently apart.

If you have any questions on the use and operation of the pelvic restraint belt, the please contact your wheelchair dealer, carer or attendant.

Notes for the patient

The pelvic restraint belt must only be fitted or adjusted by an approved SunriseMedical dealer / agent. The pelvic restraint belt must be checked on a daily basis to ensure it is adjusted correctly and free from any obstruction or adverse wear. Sunrise Medical does not recommend the transportation of a person in a vehicle with the pelvic restraint belt as a restraint system. See the transit brochure from Sunrise Medical for further information on transportation.

Maintenance:

Check pelvic restraint belt, and securing components, at regular intervals for any sign of frays or damage. Depending on usage, it may need to be replaced.

NOTE:

The pelvic restraint belt should be adjusted to suit the end user as detailed above Sunrise Medical also recommends that the length and fit of the belt is checked on a regular basis to reduce the risk of the end user inadvertently re-adjusting the belt to an excessive length.

Options - Stabilising bar

Folding stabilising bar

This bar is used to stabilise the backrest. To be able to fold the wheelchair, the release lever must be pushed inwards (Fig. 32) or released and the stabilising must be flipped downwards. When unfolding the wheelchair, please make sure that the stabilising bar is locked into position.

Fig. 32

Options - Transit wheels

Transit wheels

Travel wheels should be used wherever your wheelchair would be too wide if the rear wheels were used. After the rear wheels have been removed with the aid of the quick-release axles, the transit wheels can immediately be used to continue riding. The transit wheels are mounted so that they are approx. 3 cm above the ground when not in use. They are thus out of the way when riding, transporting, or when tipping to pass over obstacles.

CAUTION!

Without rear wheels, your wheelchair has no wheel locks.

Fig. 33

Fig. 34

Security checking and maintenance

As the user, you are the first person to notice any possible defects. We therefore recommend that before each use, you check the items in the following list:

- correct tyre pressure
- that the brakes work correctly
- that all removable parts are secure (e.g. armrest, footrest, quick-release axle...)
- no visible damage e.g. on the frame, backrest, seat sling and back sling, wheels, footplate....)

If there is any damage or malfunction, please contact your authorised dealer.

Tyres and Mounting

Tyres and Mounting

Always make sure that you maintain the correct tyre pressure (min. 3.5 bar) as this can have an effect on wheelchair performance. If the tyre pressure is too low, rolling resistance will increase, requiring more effort to move the chair forward. Low tyre pressure also has a negative impact on manoeuvrability. If the tyre pressure is too high, the tyre could burst.

The correct tyre pressure for a given tyre is printed on the surface of the tyre (min. 3.5 bar).

Tyres can be mounted the same way as an ordinary bicycle tyre. Before installing a new inner tube, you should always make sure that the base of the rim and interior of the tyre are free of foreign objects.

Check the pressure after mounting or repairing a tyre. It is critical to your safety and to the wheelchair's performance that regulation air pressure be maintained and that tyres are in good condition.

Maintenance and care

Maintenance

- Check the tyre pressure every 4 weeks. Check all tyres for wear and damage.
- Check the brakes approximately every 4 weeks to make sure that they are working properly and are easy to use.
- Change tyres as you would an ordinary bicycle tyre.
- All of the joints that are critical to using your wheelchair safely are self-locking nuts. Please check every three months to make sure that all bolts are secure (see the section on torque). Safety nuts should only be used once and should be replaced after use.
- Please use only mild household cleansers when your wheelchair is dirty. Use only soap and water when cleaning the seat sling.
- If your wheelchair should ever get wet, please dry it afterwards.
- A small amount of sewing-machine oil should be applied to quick-release axles approximately every 8 weeks. Depending on the frequency and type of use, we recommend taking your wheelchair to your authorised dealer every 6 months to have it maintained by trained personnel.

CAUTION:

Sand and sea water (or salt in the winter) can damage the bearings of the front and rear wheels. Clean the wheelchair thoroughly after exposure.

The following parts can be removed and sent back to the manufacturer/dealer for repair:

- Rear wheels:
- Armrest
- Footrest holder
- Anti-Tip Tubes

These components are available as spare parts. For further information, please see the spare parts catalogue.

Hygiene measures when being re-used:

Prior to the wheelchair being re-used, it must be carefully prepared. All surfaces which come into contact with the user must be treated with a disinfection spray.

To do this, you must use a disinfectant from the DGHM list, e.g. Antifect Liquid (Schülke & Mayr) for rapid alcohol-based disinfection for medical products and medical devices, which must be disinfected quickly.

Please take into account the manufacturer's instructions for the disinfectant you are using.

In general, a complete disinfection cannot be guaranteed on seams. We therefore recommend that you dispose of seat and back slings to avoid microbacterial contamination with active agents according to § 6 infection protection law.

Trouble Shooting

Wheelchair pulls to one side

- Check tyre pressure
- Check to make sure wheel turns easily (bearings, axle)
- Check angle of castors
- Check to make sure both castors are making proper contact with the ground

Castors begin to wobble

- Check angle of castors
- Check to make sure all bolts are secure; tighten if necessary (see the page on torque)
- Check to make sure both castors are making proper contact with the ground

Wheelchair does not fit in the support

- The chair is still new, or the seat or back slings are still stiff, this will ease through frequent use

Wheelchair is difficult to fold up

- Adjustable backrest sling is too stiff, loosen slightly.

Wheelchair squeaks and rattles

- Check to make sure all bolts are secure; tighten if necessary (see the page on torque)
- Apply small amount of lubrication to spots where movable parts

come in contact with one another

Wheelchair wobbles

- Check angle of castors
- Check tyre pressure
- Check to see if rear wheels are adjusted differently

Disposal / Recycling of materials

If the wheelchair has been made available to you free of charge, then it does not belong to you. If it is no longer required, then follow the instructions to return it as given by the organisation that made the wheelchair available to you.

In the following section, there is a description of the materials used on the wheelchair, in view of the disposal or recycling of the wheelchair and its packaging.

There may also be special local regulations in force with regard to disposal or recycling, these must be taken into account when disposing of your wheelchair. (This can include the cleaning or decontamination of the wheelchair prior to disposal).

Aluminium: Castor forks, wheels, chair sideguards, armrest frame, footrest, push handle

Steel: securement points, quick-release axle

Plastic: Handles, tube stoppers, castors, footplates, armpads and tyres

Packaging: Plastic bags made of soft polyethylene, cardboard

Upholstery: Polyester fabric with PVC coating and highly-flammable foam.

Disposal or recycling should be carried out by a disposal company or at a public disposal point. You can also return your wheelchair to your dealer for disposal.


Nameplate

Nameplate

The nameplate is located on either the cross-tube assembly or the transverse frame tube, as well as on a label in the owner's manual. The nameplate indicates the exact model designation and other technical specifications. Please provide the following pieces of information whenever you have to order replacement parts or to file a claim:

- Serial number
- Order number
- Month/Year

EASY life
072000

max 10°

40 cm

42 cm

Max
125 kg

Read the user manual

08.09.2010

Build date

SN EL0911010

Serial number

Product name, SKU number

Maximum safe slope with anti tips, depends on wheel chair setting, posture and physical capabilities of the user

Seat width

Seat depth

Maximum user weight

Read the user manual

Build date

Serial number

Technical Specifications

Technical Specifications

Overall width:

With standard wheels including handrims, close mount:

- in combination with desk sideguard: SW + 20 cm
- in combination with aluminium sideguard: SW + 18 cm
- in combination with plastic clothes protector: SW + 19 cm
- With extra-narrow handrim fitting, overall width is reduced by 2 cm

Folded dimensions:

- with rear wheels approx. 30 cm
- without rear wheels approx. 25 cm

Weight in kg:

- Transportation (without footrest, wheels, sideguard) 7 kg
- Footrest (unit) 0.8 kg
- Sideguard (unit) 1.3 kg
- Wheels (unit) 2.2 kg

Maximum user weight:

- Easy Life up to a load of 125 kg

Fig. 35

Fix Front



Fig. 36

Swing Away

Technical Specifications - Possible seat-height adjustments

According EN12183: 2009 all there listed parts (seat - back upholstery, arm pad, side panel...) are fire resistant according EN1021 Part 1/2
Maximum user weight (weight of the test dummy): 125 kg

	Min.	Max.		Min.	Max.
Overall length with footrest	935 mm	1075 mm	Seat surface angle	0°	12°
Overall width	540 mm	720 mm	Effective seat depth	360 mm	500 mm
Folded length	935 mm	1075 mm	Effective seat width	360 mm	500 mm
Folded width	340 mm	340 mm	Seat height at the front edge	430 mm	520 mm
Folded height	970 mm	970 mm	Backrest angle	-9°	9°
Total weight	11 kg	18 kg	Back height	250 mm	475 mm
Weight of the heaviest individual part	-	2.2 kg	Distance from the footrest to the seat	340 mm	470 mm
Static stability - downhill *	10°	10°	Angle from leg to seat	110°	120°
Static stability - uphill (with anti-tip tube) *	10°	10°	Distance from the armrest to the seat	230 mm	320 mm
Static stability - sideways *	10°	10°	Front position of the armrests	250 mm	360 mm
			Handrim diameter	535 mm	540 mm
Minimum turning radius	830 mm	860 mm	Horizontal axle position	230 mm	830 mm

* Maximal safe slope with anti-tips depends on wheelchair settings, posture and physical capabilities of the user

Possible seat-height adjustments Easy Life

CAUTION! measured without seat cushions





Castor - Type	Fork - Type	Seat height front in cm	only 22" Rear wheel	Seat height rear in cm 22" and 24" Rear wheel	only 24" Rear wheel
4" solid rubber	Short (4 hole) Short (4 hole) Short (4 hole) Long (5 hole) Long (5 hole)	43 45 46 47 48	40 41 40 41 40 41 40 41 40 41	43 43 44 43 44 46 43 44 46 47 43 44 46 47	48
5" solid rubber 5" soft wheel	Short (4 hole) Short (4 hole) Short (4 hole) Long (5 hole) Long (5 hole) Long (5 hole)	44 45 46 47 48 49	40 41 40 41 40 41 40 41 40 41 40 41	43 44 43 44 43 44 46 43 44 46 47 43 44 46 47 43 44 46 47	48 48 49
6" solid rubber	Long (5 hole) Long (5 hole) Long (5 hole) Long (5 hole)	47 48 49 50	40 41 40 41 40 41 - 41	43 44 46 47 43 44 46 47 43 44 46 47 43 44 46 47	48 48 49 48 49
6" soft wheel	Long (5 hole) Long (5 hole) Long (5 hole) Long (5 hole)	46 48 49 50	40 41 40 41 40 41 - 41	43 44 46 43 44 46 47 43 44 46 47 43 44 46 47	48 48 49 48 49
7" solid rubber	Long (5 hole) Long (5 hole) Long (5 hole) Long (5 hole)	48 49 50 52	40 41 40 41 41 - -	43 44 46 47 43 44 46 47 43 44 46 47 43 44 46 47	48 48 49 48 49 48 49 51

The wheelchair conforms to the following standards:





- a) requirements and test methods for static, impact and fatigue strengths (ISO 7176-8)
b) Flammability resistance of upholstered parts in accordance with ISO 7176-16 (EN 1021-1/2)

Yes
Yes

Torque

Torque for M6 screws is 7 Nm, unless otherwise specified.

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Sunrise Medical AB
Sunrise Medical B.V.
V. Guldman A/S - G2